

# Accutome AccuPen®

## Handheld Applanation Tonometer

Accutome AccuPen 24-3000

Busy facilities deserve an instrument that will come through in the clutch. The AccuPen's combination of accuracy and portability will make it an invaluable instrument to your entire staff.

### • Accurate

Gravity Offset Technology provides precise IOP measurements with less calibration compared to other handheld tonometers.

Micro strain gauge technology combined with proprietary algorithm provides reproducible results.

### • Easy-to-use

Ergonomic design combined with sighting lines allow for easy visualization of the cornea, facilitating both centralization and perpendicularity.

Contoured molding with large, visible LCD allows the AccuPen to be used ambidextrously and fit comfortably in any size hand.

### • Ready-to-use

Eliminates the constant need for "up, down and around" calibration.

Long lasting lithium battery takes up to 15,000 readings.

### • Portability

Weighs only 3.0 oz.

Comes with custom, heavy duty traveling case and lanyard for safe transportation.



### • Versatility

Takes fast measurements with patient in any position.

Frees up valuable office space.

### • Smarter Software

Stores and averages nine measurements.

Built in IOP Correction Calculator adjusts IOP based on corneal thickness.

### Specifications:

Accuracy\* +/- 2mmHg

Weight - 3.0 oz

Dimensions - 7 1/4" X 1 1/4"

Transducer Dimensions - 1.0 mm diameter

Measurement Range - 5 - 60mmHg

Power Source - 3.6 V Lithium Battery

\*Compared to Goldmann Tonometer



Accutome Repair Center | Toll Free: 800-979-2020

24 hour turnaround for most brands of Handheld instruments and Ultrasound Equipment

Optos, Inc.  
67 Forest Street  
Marlborough, MA 01752  
USA

For sales questions and orders: 800-854-3039

©2011 Distributed by Optos. All rights reserved. Optos® plc. P/N 433001/1



optos.com